

## Recommended Starting Material

The size of recommended starting material to use with determined lysis volume are listed here.

**Table:** Recommended starting material and Lysis Buffer amount

Sample	Size of Starting Material	Lysis Buffer Amount
Under 18 weeks	5-6 ml of amniotic fluid	250 µl
Above 18 weeks	2-3 ml of amniotic fluid	

## Protocol

### *Isolation of Genomic DNA (based on silica technology)*

#### Sample type:

- Amniotic fluid

#### Some tips to know:

- All centrifugation steps are carried out at room temperature (15–25°C).
- If RNase treatment is desired, Prime-RNase can be ordered separately from ROJETechnologies, Cat NO. EB983013.
- Phosphate buffered saline (PBS) can be prepared as written in appendix 4 in DNJia AmnioPure Kit or can be ordered separately from ROJETechnologies, Cat NO. BU983026.
- Do not forget to add the appropriate amount of Ethanol (96–100%) to BW1 and BW2 buffers (BWB1 and BWB2) as indicated on the bottle, before using for the first time.
- If AL Buffer (ALB) forms precipitate, please warm it to 56°C until the precipitate has fully dissolved.

## Process

1. Centrifuge 1.5 ml of the amniotic fluid for 10 min at 5000 RPM in a clean microcentrifuge tube, discard the supernatant and repeat this step until you get pellet from whole sample mentioned in the table (based on the gestation week).
2. Add 20  $\mu$ l RJ-Protease to the pellet. Add 250  $\mu$ l PBS to cell pellets to the tube. Then add 250  $\mu$ l ALB. Pulse vortex for 15 s and incubate at 56 °C for 25 min.

**Optional:** Pulse-vortexing for 3 times during incubation may increase the DNA yield.

**Optional:** This protocol is specialized for DNA isolation. For usual and routine applications RNase treatment isn't needed. If RNA-free genomic DNA is required, add 10  $\mu$ l Prime-RNase A (Cat NO. EB983013), mix by vortexing, and incubate for 5 min at room temperature before going to step 2.

3. Add 250  $\mu$ l Ethanol (96-100%, molecular biology grade) to the lysate, mix by pulse vortexing for 15 s, then centrifuge briefly.
4. Gently, pipette the mixture to a spin column placed in a 2 ml collection tube (supplied in the kit box). Centrifuge at 8000 rpm for 1 min. Discard flow-through.

**Note:** If the lysate did not pass the column, centrifuge at full speed until it passes through the column.

5. Add 600  $\mu$ l BWB1, incubate for 2 min at room temperature (15-25°C) and centrifuge for 1 min at 8000 rpm. Then discard the flow-through.
6. Add 600  $\mu$ l BWB2, incubate for 2 min at room temperature (15-25°C) and centrifuge for 3 min at 14000 rpm. Discard the flow-through and centrifuge for 1 min at 14000. Then discard both the flow-through and the collection tube. Place the spin column in a new clean 1.5 ml microcentrifuge tube (not provided).
7. Pipette 25-30  $\mu$ l RRB directly onto spin column. Incubate at 37°C for 1-5 min. Centrifuge it at 8000 rpm for 1 min.

**Note:** If the expected DNA yield is more than the yield from pervious step, put the spin column on a new microtube and add another 25-30  $\mu$ l RRB, incubate for 5 min at 37°C. Then, centrifuge for 1 min at 8000 rpm.