

Hashin, 6X

DNA Loading Dye, Blue

Store at 15-25 °C

- 2ml: Cat. #LD983003
- 5ml: Cat. #LD983004

Description

Hashin, 6X is a ready to-use buffer for loading DNA samples. This blue loading dye is ideal for amplicons larger than 200 bp on 0.5 to 1.5 percent agarose gel and larger than 100 bp on 2 to 3 percent agarose gel.

Note Add 1/6 volume of 6X DNA Loading Dye to your DNA sample.

Quality Control

Hashin, 6X is tested against predetermined experiments on a lot-to-lot basis according to ROJETechnologies ISO-certified quality management system, to ensure consistent product quality. For your information, the results of all experiments are accessible by addressing REF and Lot number on web at www.rojetechnologies.com.

Technical Assistance

ROJETechnologies guarantees your complete satisfaction. ROJE technical support team composed of highly trained experienced scientists, who are able to troubleshoot most problems you face. Our technical support team can offer expert advice which may help you select suitable product.

Contact our technical support at any time by selecting one of these ways:

- Through our telephone and fax number; +98 353 7237122, +98 353 7302468.
- You can submit your question directly to ROJE Technical Support Team from our website (www.rojetechnologies.com).
- Or send your questions to this email address, technicalsupport@rojetechnologies.com .

For Research Use Only. Not for use in diagnostic procedures.

By ROJE, Edition 01/2022

www.rojetechnologies.com