

RT-PCR Grade Water

Description

ROJE RT-PCR Grade Water is an ultrapure nuclease-free water perfectly ideal to be used in RT-PCR reactions. This product is tested for RT-PCR applications in a lot-to-lot basis. NOT DEPC treated.

Features

Specific parameters of RT-PCR Grade Water are listed here in Table 1.

Table 1: RT-PCR Grade Water specifications

Parameters	
Formulation	H ₂ O
Appearance	Colorless, one uniform clear phase
Toxicity	Standard handling
Storage	+15°C to +25°C
Do not freeze	Ok to freeze
Endotoxins	<0.25 IU/ml
TOC content	<500 ppb
RNase/ DNase	0
Sterility	+

Ordering Information

Category	Product name	Cat NO.	Size
Water, Buffers and Solutions	RT-PCR Grade Water	WA983020	1.5ml

Technical Assistance

ROJETechnologies guarantees your complete satisfaction. ROJE technical support team composed of highly trained experienced scientists, who are able to troubleshoot most problems you face. Our technical support team can offer expert advice which may help you select suitable product.

Contact our technical support at any time by selecting one of these ways:

- Through our telephone and fax number available at ROJETechnologies website.
- You can submit your question directly to ROJE Technical Support Team from our website (www.rojetechnologies.com).
- Or send your questions to this email address, Technicalsupport@Rojetechnologies.com.

For Research Use Only. Not for use in diagnostic procedures.

www.rojetechnologies.com

By ROJE

Edition: 01/2022